

Making improved predictions of Non-invasive Continuous Arterial Blood Pressure through leveraging EHR medication data

SIMON LEE¹, Ákos Rudas², & Jeffery N. Chiang²

1 BIG Summer Program, Institute for Quantitative and Computational Biosciences, UCLA

2 Department of Computational Medicine, University of California, Los Angeles, CA, USA

Background



Non-invasive non-continuous blood pressure measurements




Invasive **continuous** blood pressure measurements



**There exists a
computational tool that
can predict ABP
non-invasively and
continuously**

ABPImputation Package



Hill, B.L., Rakocz, N., Rudas, Á. et al. Imputation of the continuous arterial line blood pressure waveform from non-invasive measurements using deep learning. *Sci Rep* 11, 15755 (2021).
<https://doi.org/10.1038/s41598-021-94913-y>


Table 2. RMSE table

Method	MIMIC	UCLA
PPG scaling	6.895 (6.876–6.914)	9.108 (9.078–9.137)
Sideris et al	13.940 (13.901–13.978)	13.111 (13.072–13.151)
1D V-Net	5.823 (5.806–5.840)	6.961 (6.937–6.985)


Table 2. Root mean square error (mean (95% CI)) for each cohort.

Predictions are pretty good! but can always be better!

How can we make predictions better




Pressor	Receptor	Initial Dose	Notes
Norepinephrine	$\alpha > \beta$	7 mcg/min	Bees knees
Epinephrine	$\beta > \alpha$	5 mcg/min	More inotropy. Consider 5mcg/min + norepi
Vasopressin	V_1	0.03 U/min	Fixed dose, no titration
Phenylephrine	α	30 mcg/min	Good for tachy hearts
Dopamine	$D \rightarrow \beta \rightarrow \alpha$	5 - 10 mcg/kg/min	Low dose for stable bradycardia
Dobutamine	β	2.5mcg/kg/min	Poor mans milrinone. Inotropy alone
Milrinone	PDEi*	0.25mcg/kg/min	Bridge to transplant. Inotropy + vasodilation
Methylene Blue	NOi, GCi*	1 - 2 mg/kg bolus	Alright Dr. Smarty-pants
Angiotensin II	ATII-R	20ng/kg/min	Alright Dr. Money-bags




Our Input Data

EKG, PPG, & Medication time series




Data Scrapping (much of my summer...)



LOGIC FLOW

1. Read `D_items.csv` to extract vasopressor drug ID's
2. Filter `INPUTEVENTS_MV.csv` for only patients who recieved vasopressors
3. Scrape the amount of drug administered
4. Calculate in milliseconds, the administration rate
5. Glob.glob the specific patient files
6. Go to `p0XXXXXX.csv.gz` (patient file) and build 5 zero columns for the vasopressor features
7. Match the drugs that were administered to the corresponding columns of the patient file
8. Find start and end time and fill in the rate
9. Write out csv with new drug columns (Most painful part)

Optimization via parallel processing



1. Start at the start and end of patient list
2. Allocate 10 gpu cores each and begin writing out csv files with new drug columns
3. Big Idea is widely inspired by TV show “Silicon Valley” - middle out compression algorithm



Complications along the way...

Monday, July 18th ▾



simonlee711 8:27 AM

I've been trying to log onto the halperin gpu server for some time and its having trouble with this issue:

```
Could not establish connection to "halperingpu01.cs.ucla.edu": Connecting with SSH
timed out.
```



Jeff Chiang 11 days ago


iget that as well...





I ended up crashing the halperin gpu server because of exceeding RAM capabilities. But finally generated my enormous dataset

Current Patient sample



- 1 Patient
 - Split into training frame and testing frames
 - (32408, 400, 15)







Bland Altman on 1D Vnet


Bland-Altman Systolic ABP - Test: 3.31 ± 7.60 Bland-Altman Diastolic ABP - Test: -3.50 ± 6.15



Bland Altman on 1D Vnet with Medication

Bland-Altman Systolic ABP - Test: 4.51 ± 6.30 Bland-Altman Diastolic ABP - Test: -3.16 ± 4.16





Q-Q normal plot

